



how it works

IdealClip is applied by simply pressing down with the fingertips and it opens up easily by pulling the protruding strip of the package

how to close









IdealClip is an airtight and reusable solution for easy reclosing any kind of flexible bags, pouches and flow-packs. Easy to apply and to open, IdealClip keeps the freshness of the product intact once the package is opened and offers and unbelievable watertight seal. It's available in different sizes, shapes and colours and can be personalised

Airtight, watertight seal

easy airtight closure

> immediate automatic opening

available in different sizes shapes and colours







STANDARD

clip suitable for reclosing any types of flexible bags and pouches



IdealClip can be inserted in-line inside the package







IdealClip packed in flow-pack can be sealed outside



STANDARD

scala 1:1





Patent	International patent
Material	Polypropylene (PP) e Polyetilene (PE)
Sizes (mm)	70 / 100 / 140 / 180
Weight (gr)	2 / 3 / 4,3 / 5,5
Standard colours	Any Pantone colours
Segment of application	Food, non-food and medical
Personalisation (embossed logo)	Only upper side (mm) 4x25 / 4x50



TWISTING

clip reinforced suitable for reclosing twisting flow-pack











TWISTING

scala 1:1





Patent	International patent
Material	Polypropylene (PP)
Sizes (mm)	35
Weight (gr)	2
Standard colours	Any Pantone colours
Segment of application	Food, non-food
Personalisation (embossed logo)	Only upper side (mm) 4x25